



Aerotel
LifeCare

SKeeper™

Personal Wristop Cellular Communicator with Emergency Response

Peace of Mind is Everywhere

Sound the alarm, locate and communicate. Anytime, anywhere. That's the key to Aerotel's SKeeper – a "peace of mind" product line designed for elderly, chronically ill, children or lone workers, as well as their relatives and caregivers.

A New Level of Personal Safety

SKeeper has taken PERS (Personal Emergency Response System) and social alarm systems, to the next level. Fitted with a distress button and built-in cellular communicator, it can activate a warning signal when needed, enabling wearers to get immediate care and help from a telecare monitoring center.



Built-in Cellular Speakerphone

Using its embedded GSM module and built-in speakerphone, SKeeper enables cellular voice calls to be made to pre-defined numbers (e.g. a relative or a family doctor) or to be received from any caller or from the remote monitoring center when in need.

Online Location and Tracking (optional service)

SKeeper can take advantage of operator provided location-based services. In the event that the wearer wanders outside a specified zone (e.g. neighbourhood or school area), the system can send an alert to the monitoring center and/or to another mobile phone.

Technical Data

Dimensions	65 x 46 x 17 mm	2.6" x 1.8" x 0.7"
Weight (battery included)	70g	2.5 oz
Operating temperature	0 ÷ 55 °C	14 ° ÷ 130 °F
Battery	3.7 VDC/630mAh Li-Ion (rechargeable)	
Standby time (typical)	Up to 100 hours	
Talking time (typical)	1 hour	
Safety call button	1	
Speed dialing buttons	3 (user programmable)	
Speaker and microphone	Internal	
Automatic call answering	Included (programmable)	
Wireless module	Cinterion MC55/MC56	
Network	Dual band GSM/GPRS (900/1800 or 850/1900)	
Antenna	Integrated	

Specifications are subject to change without notice.

Aerotel Medical Systems (1998) Ltd.

5 Hazoref St. Holon 58856, Israel

Tel +972 3 559 3222

Fax +972 3 559 6111

info@aerotel.com



www.aerotel.com